TIPSAN TIBBİ ALETLER SAN. VE TİC. A.Ş.  
Kemalpaşa Cad.7404/1 Sk.No:3 Pınarbaşı- İzmir- Turkey  
Tel: +90 232 479 5654 www.tipsan.com.tr

**SURGICAL INSTRUMENTS**  
**TPS-IFU-03 SURGICAL INSTRUMENTS**

Rev.00, Issue Date:10.10.2022

The following languages are included in this packet:

English (EN)      Turkish (TR)

For additional information and translations please contact the manufacturer or local distributor.

*Attention Operating Surgeon*  
**IMPORTANT MEDICAL INFORMATION**  
**SURGICAL INSTRUMENTS**

- A.DEFINITIONS
- B.DESCRPTION
- C.INTENDED PURPOSE & INDICATIONS
- D.CONTRAINDICATIONS
- E.GENERAL PRECAUTIONS &WARNINGS
- F.ADVERSE EFFECTS
- G. LIMITATIONS AND RESCTRCTIONS OF REPROCESSING
- H. CLEANING /DISINFECTION
- I. INSPECTION, MAITENANCE AND TESTING
- K. PACKAGING
- L. STERILIZATION
- M. STORAGE
- N.SUBSTANCES CONTAINED BY THE DEVICE
- O.REPORTING OF SERIOUS INCIDENTS

**A.DEFINITIONS:**

Symbols and abbreviations may be used on the package label. The following table provides the definition of these symbols and abbreviations.

Table 1. Definitions of Symbols and Abbreviations

Symbol	Definition	Symbol	Definition
	Batch code		Manufacturer
	Catalog number		Unique Device Identifier Number
	Caution, consult accompanying documents		Medical Device
	Consult operating instructions		Non-sterile
	Date of manufacture		

**B.DESCRPTION:**

The Surgical instruments Supplied by Tipsan are reusable orthopaedic surgical instruments designed solely for use in fitting of Tipsan orthopaedic implants.

This instruction manual pertains to all hip, knee, bone fixation, and extremity reusable medical devices manufactured and distributed by Tipsan

Surgical instruments are supplied non-sterile and must be cleaned and sterilized before use. After use, these instruments must be, properly decontaminated, cleaned, and stored. The following information outlines the proper steps for reprocessing Tipsan surgical instruments to help assure their long life. The Surgical instruments including drills are not intended to be connected to an active device.

These instructions are recommended for the care, cleaning, maintenance and sterilization of reusable Tipsan Orthopaedic surgical instruments. This document is intended to assist health care personnel in safe handling practices, effective reprocessing and maintenance of Tipsan reusable devices. It provides information complementary to the instructions for use in fulfilment of EN ISO 17664, ISO 16061 and the Regulation (EU) 2017/745 of the European Parliament and of the Council, Annex 1, section 23.4(n).

**Materials:**

The Surgical Instruments are manufactured from a variety of materials which include Stainless Steel, titanium alloy, Aluminium and plastic materials all of which conform to ASTM or ISO standards, or internal standards.

**C.INTENDED PURPOSE & INDICATIONS****Intended Purpose:**

Tipsan Surgical Instruments are designed to use to help surgeon to place orthopaedic implants to adult patient properly and facilitate bone preparation.

**Indications:**

the Surgical Instruments are indicated for use in orthopaedic Knee, Hip replacement and Bone Fixation implant surgery.

**Intended Clinical Performance:**

Intended Clinical performance of device is:

- place implants to patient properly.

**Intended Patient populations:**

The Device is designed to be used for Adult Patients . No limitations in regard to gender or aetiology

**Intended user :**

The Device is designated to be used only by Health care professional – Orthopaedic Surgeon,

**D.CONTRAINDICATIONS**

The Instruments should not be used for anything other than their intended use

**E.GENERAL PRECAUTIONS & WARNINGS**

1- Under no circumstances may Tipsan Instruments be used to fit components from other manufacturers, if Tipsan instruments are used in such

cases the company shall not be liable for the operation of the resulting unit. Likewise, Tipsan implants may only be implanted with the appropriate instruments supplied by Tipsan. Tipsan shall not be liable for implants implanted using other manufacturer's instruments

- 2- As a result of mechanical features required, most of the instruments are made of non implantable material. In the event an instrument breaks, no fragment must remain in the patient as this could cause post-operative complications such as allergies, infections, or complications of a biological nature associated with the release of metal components, possibly requiring further intervention
- 3- The life of the instruments depends on the number of times they are used as well as the precautions taken in handling, cleaning and storage. Great care must be taken of the instruments to ensure that they remain in good working order
- 4- Instruments should be examined for wear or damage by Surgeon and staff in operating centers prior surgery, the examination shall include a visual and functional inspection of the working surfaces, articulation points, and springs. It should also include verifying all welded connections, that all components are present, and cleanliness of the orifices and cavities, as well as the absence on any cracks, distortion, impact, corrosion other change
- 5- Any complications or other effects that may occur for reasons such as incorrect indication or surgical technique, unsuitable choice of material or treatment, inappropriate use or handling of the instruments etc., fall under the responsibility of the operating surgeon and can not be blamed on the manufacturer
- 6- The Manufacturer does not accept liability in the case of non-compliance with the requirement of this user instruction.
- 7- Instruments that are contaminated, unsterile or damaged or that have been improperly handled or instruments or their components that have been altered without authorization must not be used under any circumstances.
- 8- If instruments are subjected to excessive loading, are damaged or are handled improperly, they may fracture, corrode, oxidize, become loose and show excessive wear or their functionality may be impaired.
- 9- Instruments may sometimes have sharp edges or pointed tips. Consequently, if plastic gloves are used there is a danger that these may become damaged. Therefore please pay attention to the risk of infection.

#### F.ADVVERSE EFFECTS

In any surgical procedure, the potential for complications exists. The risks and complications with these surgical instruments include:

- 1- Infection, if the instruments are not properly cleaned and sterilized
- 2- Breakage of instruments
- 3- Cutting the gloves or skin of the surgical staff
- 4- Crack, fracture of the bone

#### G. LIMITATIONS AND RESCTRCTIONS OF REPROCESSING

Surgical instruments are designed for their durability and ability for reuse. Tipsan's reusable instruments are typically manufactured from stainless steel and polymers, which permits a long life when handled and maintained properly. Repeated processing has minimal effect on these instruments. End of functional life is normally determined by wear and damage due to use. Tipsan does not define a maximum number of reuse/reprocessing cycles for the reusable instruments, The life span of these instruments depends on how they are cared for during use, cleaning and handling

Broken/damaged instruments should be discarded in accordance with all local and national regulations for infection prevention and control.

#### H. CLEANING

The Cleaning Processes (automatic ) presented in this instruction have been validated. Other methods of cleaning may be suitable but must be validated by the user of the device. The manual method is effective for all devices may be used when an automated option is not available

**Note:** To ensure proper processing, please follow appropriate healthcare facility standards.

<b>Initial treatment at point of use</b>	<b>Clean instruments as soon as possible after use.</b> Do not allow blood or debris to dry on the instruments. If cleaning must be delayed, place groups of instruments in a covered container with cold water or an appropriate detergent or Alkaline solution to delay drying. Clean all instruments whether or not they were used or inadvertently contacted with blood or saline solution.
<b>Warnings</b>	<ul style="list-style-type: none"> <li>• When handling sharp instruments use extreme caution to avoid injury: consult with an infection control practitioner to develop and verify safety procedures appropriate for all levels of direct instrument contact.</li> <li>• Instruments must not come into contact with substances containing chlorine or fluorine. Instruments that are made wholly of plastic must not come into contact with strong acids or organic or ammonia-containing solvents, aromatic and/or halogen hydrocarbons or oxidising chemicals. Aluminium and materials containing aluminium must never come into contact with substances containing mercury. Even the smallest traces mercury can lead to considerable corrosion. Instruments made from materials containing aluminium must be wiped and cleaned with, or placed in detergents and disinfections that have a pH value between 4.5 and 8.5. with higher or lower pH values the protective neutral coating of materials containing aluminium dissolves, which leads to corrosion. Detergents containing caustic may cause the surface corrosions and discoloration on anodized aluminium Instrument set therefore Detergent containing caustic should not be used during cleaning process for those parts.</li> </ul>
<b>Preparation for Cleaning</b>	<ul style="list-style-type: none"> <li>•The cleaning process must be conducted so that all parts of the surgical instrument are exposed as permitted by instrument design. The cleaning process should include an individual properly gowned with appropriate glove and personal protective equipment.</li> <li>• This may require opening all hinged items or the disassembly of those items with multiple or removable parts.</li> <li>• Those items with mating surfaces, i.e. ratchets, hinges, serrations, lumens, blind holes, etc. must be carefully cleaned to remove all visible debris from the items.</li> <li>• Additional assembly/disassembly instructions may be found in the product specific surgical technique.</li> </ul>

<b>Automated Cleaning</b>	<p><b>Required Equipment:</b> Automatic Washer machine: an Automatic washer machine with approved efficacy should be used ( CE mark and validated according to ISO 15883), properly installed, qualified and regularly subjected to maintenance and testing</p> <p>Cleaning Accessories: Brushes and/or Pipe Cleaners, Syringes, Gloves, Absorbent Disposable Cloth</p> <p>Note:use Cleaning detergent intended for use in automated cleaning, Do not exceed the concentration and temperature recommended by the detergent manufacturer</p>																											
	<p>Manual Precleaning is not required if the device does not have dried-on visible debris. If the device does have dried on visible debris, follow the manual cleaning steps prior to placing in the device in the automatic washer for cleaning.</p> <ol style="list-style-type: none"> <li><b>Separate</b> any mated instruments before cleaning. For any instruments with moving pieces, move the pieces throughout their range of motion during cleaning to clean moving pieces in all positions.</li> <li><b>Rinse</b> with cold tap water (&lt; 45°C) to remove visible contamination. While rinsing, scrub thoroughly with a soft brush and/or pipe cleaner and repeatedly flush lumens and blind holes with a syringe.</li> <li><b>Transfer</b> to washer for processing. See table below for cycle parameters</li> </ol> <p>Washer Parameters</p> <table border="1"> <thead> <tr> <th>Phase</th> <th>Recirculation Time (Minutes)</th> <th>Temperature</th> <th>Detergent/Water Type</th> </tr> </thead> <tbody> <tr> <td>Pre-washing</td> <td>2</td> <td>30 °C</td> <td>Deionized water</td> </tr> <tr> <td>Cleaning</td> <td>7</td> <td>55°C</td> <td>alkaline-based Detergent, (Dosage:120 ml/l) For Example: METALCLEAN ; Ecolab Ltd. - Deionized Water</td> </tr> <tr> <td>Rinsing</td> <td>1</td> <td>55°C</td> <td>Deionized Water</td> </tr> <tr> <td>Final Rinsing</td> <td>5</td> <td>55°C</td> <td>Deionized Water</td> </tr> <tr> <td>Disinfection</td> <td>5</td> <td>90°C</td> <td>Deionized Water</td> </tr> <tr> <td>Drying</td> <td>25</td> <td>80°C</td> <td>N/A</td> </tr> </tbody> </table> <ol style="list-style-type: none"> <li><b>Visually inspect</b> for cleanliness. All visible surfaces, internal and external, should be visually inspected. If necessary, re-clean until it is visibly clean.</li> </ol> <p>Note: The process above is validated; however, variations in process parameters or equipment may compromise the sterility assurance level.</p>	Phase	Recirculation Time (Minutes)	Temperature	Detergent/Water Type	Pre-washing	2	30 °C	Deionized water	Cleaning	7	55°C	alkaline-based Detergent, (Dosage:120 ml/l) For Example: METALCLEAN ; Ecolab Ltd. - Deionized Water	Rinsing	1	55°C	Deionized Water	Final Rinsing	5	55°C	Deionized Water	Disinfection	5	90°C	Deionized Water	Drying	25	80°C
Phase	Recirculation Time (Minutes)	Temperature	Detergent/Water Type																									
Pre-washing	2	30 °C	Deionized water																									
Cleaning	7	55°C	alkaline-based Detergent, (Dosage:120 ml/l) For Example: METALCLEAN ; Ecolab Ltd. - Deionized Water																									
Rinsing	1	55°C	Deionized Water																									
Final Rinsing	5	55°C	Deionized Water																									
Disinfection	5	90°C	Deionized Water																									
Drying	25	80°C	N/A																									

## I. INSPECTION, MAINTENANCE AND TESTING

Surgical instruments and instrument cases are susceptible to damage from prolonged use, and through misuse or rough handling. Care must be taken to avoid compromising their exacting performance. To minimize damage, the following should be done:

- Inspect the instrument case and instruments for damage when received and after each use and cleaning. Incompletely cleaned instruments should be re-cleaned, and those that need repair should be returned to Tipsan.
- After cleaning, the disassembled instruments should be reassembled and placed in their proper locations in the instrument cases where appropriate. Locking mechanisms should fasten securely.
- Only use an instrument for its intended purpose.
- For devices with hinged or moving components; a biocompatible, surgical-grade lubricant intended for heat sterilized medical instruments should be used.

Tipsan does not accept responsibility or liability of this instrument nor any of the component parts upon which repairs and/or modifications have been made or attempted except as performed by Tipsan

The surgical instruments are not subjected to periodic maintenance and calibration process.

## K. PACKAGING

Tipsan instrument cases are intended to protect instrumentation during shipping. Health care personnel bear the ultimate responsibility for ensuring that any packaging method or material, including a reusable rigid container system, is suitable for use in sterilization processing and sterility maintenance in a particular health care facility. Testing should be conducted in the health care facility to assure that conditions essential to sterilization can be achieved. Tipsan does not accept responsibility or liability arising from a lack of cleanliness or sterility of any medical devices supplied by Tipsan that should have been cleaned and sterilized by the end user.

## L. STERILIZATION

Tipsan instruments manufactured of stainless steel may be steam sterilized with no detrimental effects. plastics can be steam sterilized. All items to be sterilized must be thoroughly cleaned and packaged appropriately for the type of sterilization. The package must permit contact of the sterilant with the item, while also serving as a barrier to microorganisms, during any storage period. Users should wear non-linting gloves, i.e. Latex or Nitrile, when handling reusable instruments, to minimize bioburden and particulates. Inspect the product packaging for tears, holes, moisture or other defects. If these concerns are present, segregate these items and reprocess them.

The following process parameters are validated by Tipsan and recommended for sterilization.

1-Double wrap the component in a non-woven medical grade wrapping material.

2-Autoclave according to the following parameters

Steam Sterilization		
Cycle Type	Parameter	Minimum Set Point
Prevacuum 273°F (134°C)	Exposure Temperature	273°F (134°C)
	Exposure Time	5 minutes

3-These recommendations have been developed and validated using specific equipment. Due to variations in environment and equipment, it must be demonstrated that these recommendations produce sterility in your environment.

#### **M. STORAGE**

All instruments must be stored in a clean, dry environment

#### **N.SUBSTANCES CONTAINED BY THE DEVICE**

the devices do not contain or incorporate:

- a medicinal substance, including a human blood or plasma derivative, or
- tissues or cells, or their derivatives, of human origin, or
- tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No 722/2012; or

#### **O.REPORTING OF SERIOUS INCIDENTS**

Surgeons/Health Facilities must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the serious incident has occurred;

TIPSAN TIBBİ ALETLER SAN. VE TİC. A.Ş.  
Kemalpaşa Cad.7404/1 Sk.No:3 Pınarbaşı- İzmir- Turkey  
Tel: +90 232 479 5654 www.tipsan.com.tr

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- (AB) 722/2012 sayılı Tüzük'te belirtildiği şekilde, hayvan kaynaklı dokular veya hücreler ya da onların türevlerini;