TIPSAN TIBBİ ALETLER SAN. VE TİC. A.Ş.  
Kemalpaşa Cad.7404/1 Sk.No:3 Pınarbaşı- Izmir- Turkey  
Tel: +90 232 479 5654 www.tipsan.com.tr

## TIPSAN BONE PLATE&SCREW SYSTEM

### TPS-IFU-05

Rev.01, Issue Date:20.05.2024

The following languages are included in this packet:

The following languages are included in this packet:

English (EN) Turkish (TR)

For additional information and translations please contact the manufacturer or local distributor.

---

*Attention Operating Surgeon*  
**IMPORTANT MEDICAL INFORMATION**  
**BONE PLATE&SCREW SYSTEM**

- A- DEFINITIONS
- B- DESCRIPTION
- C- INTENDED PURPOSE & INDICATIONS
- D- CONTRAINDICATIONS
- E- GENERAL PRECAUTIONS & WARNINGS
- F- ADVERSE EFFECTS
- G- HANDLING AND STERILIZATION
- H- STORAGE CONDITIONS
- I- MAGNATIC RESONANCE IMAGING (MRI) SAFETY
- J- SUBSTANCES CONTAINED BY THE DEVICE
- K - REMOVEL AND DISPOSAL
- L- IMPLANT CARD
- M- REPORTING OF SERIOUS INCIDENTS
- N- OTHER INFORMATION

**A.DEFINITIONS:**

Symbols and abbreviations may be used on the package label. The following table provides the definition of these symbols and abbreviations.

Table 1. Definitions of Symbols and Abbreviations

Symbol	Definition	Symbol	Definition
	Batch number		Date of manufacture
	Catalog number		Manufacturer
	Do not re-use		Unique Device Identifier Number
	Caution, consult accompanying documents		Medical Device
	Consult operating instructions		Non-sterile

**B.DESCRPTION:**

**Tipsan** Fracture fixation devices (Plate& Screw System) are used only as an aid to healing; they are not a substitute for normal intact tissue or bone. The anatomy of human bones presents limitations with respect to the size or thickness of bone screws or plates and thus the strength of implants is limited. Fracture fixation devices are available in many styles and sizes and are made from various types of metal used in surgical implants. The component material is provided on the implants and labels. Do not mix components from different manufacturer. All implantable devices are designed for single use only. This Fracture fixation devices are not approved for fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

Through the advancement of surgical fusion hardware, the surgeon has been provided a means of correcting deformity and reducing pain for many patients. While the implants used are largely successful in attaining these goals, it must be recognized that they are manufactured from metal, and that no implant can be expected to withstand the activity levels and loads as would normal, healthy bone after fusion occurs. The surgeon must evaluate each situation individually based on the patient’s clinical presentation in making any decisions regarding implant selection.

The implants are single use only devices.

**Materials:**

The implants are manufactured from a variety of materials which include Unalloyed Titanium, titanium alloy, Stainless Steel all of which conform to ASTM or ISO standards, or internal standards. The component material is provided on the package label

Table-2 : Component List & Materials & Application Method.

Components	Materials
Plates	ISO 5832-3 Titanium Alloy; ISO 5832-2 Unalloyed Titanium, ISO 5832-1 Stainless Steel
Screws	ISO 5832-3 Titanium Alloy, 5832-1 Stainless Steel

Implantation of the devices will expose patients to the following materials:

Material Description	Elements in Material	% Composition
Titanium alloy - Ti-6Al-4V ( ISO 5832-3 )	Aluminium	5,5 to 6,75
	Vanadium	3,5 to 4,5
	Iron	0,3 max.
	Oxygen	0,2 max.
	Carbon	0,08 max.
	Nitrogen	0,05 max.
	Hydrogen	0,015 max.
	Titanium	Balance

Material Description	Elements in Material	% Composition
Unalloyed Titanium ( ISO 5832-2)	Nitrogen	0,03 max.
	Carbon	0,08 max.
	Hydrogen	0,0125 max.
	Iron	0,30 max.
	Oxygen	0,25 max.
	Titanium	Balance

Material Description	Elements in Material	% Composition
Stainless Steel (ISO 5832-1)	Carbon	0,030 max.
	Silicon	0,75 max.
	Manganese	2,0 max.
	Phosphorus	0,025 max.
	Sulfur	0,010 max.
	Nitrogen	0,10 max.
	Chromium	17,0 to 19,0
	Molybdenum	2,25 to 3,00
	Nickel	13,0 to 15,0
	Copper	0,50 max.
	Cobalt	0,10 max.
	Iron	Balance

### C.INTENDED PURPOSE & INDICATIONS

#### Intended Purpose:

Plate & Screw System is intended for temporary fixation, correction or stabilization of fresh bone fractures at upper and lower extremity of adult patients.

#### Indications:

-fresh bone fractures at upper and lower extremity of adult patients.

#### Intended Clinical Performance:

Intended Clinical performance of device is:

- Pain reduction;
- Bone Fracture Healing

However, these performances are dependent on several parameters regardless of the device concerned:

- The patient (weight, activity level, expectations, age, etc.);
- The patient's bone quality
- The implantation quality.

#### Intended Patient populations:

The Device is designed to be used for Adult Patients. No limitations in regard to gender or aetiology

#### Intended user :

The Device is designated to be used only by Health care professional – Orthopaedic Surgeon,

### D.CONTRAINDICATIONS

1. Infection
2. Physiologically or psychologically inadequate patient
3. Inadequate skin, bone, or neurovascular status
4. Irreparable tendon system
5. Possibility for conservative treatment
6. Patients with high levels of activity

### E.GENERAL PRECAUTIONS & WARNINGS

-The manufacturer is not liable for complications that might occur for reasons such as incorrect indications or surgical technique, unsuitable

choice or handling of material and/or surgical instruments.

-The correct selection and sizing of the implant is extremely important. Selection of the proper size, shape, and design of the implant increases the potential for success. The implants require careful seating and adequate bone support.

-In selecting patients for surgery, the following factors can be critical to the eventual success of the procedure:

1. Patient's occupation or activity. If the patient is involved in an occupation or activity which includes substantial lifting or muscle strain, the resultant forces can cause failure of the fixation, the device, or both. The implant will not restore function to the level expected with normal healthy bone, and the patient should not have unrealistic functional expectations.
2. Condition of senility, mental illness, or alcoholism. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made prior to material selection or implantation.

-Following the instructions for use provided in product literature can minimize the potential for complications or adverse reactions with any implant.

-It is the responsibility of each surgeon using implants to consider the clinical and medical status of each patient and to be knowledgeable about all aspects of implant procedure and the potential complications that may occur. The benefits derived from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are common. The patient's mental status must also be considered. Willingness and/or ability to follow post-operative instructions may also impact the surgical outcome. Surgeons must balance many considerations to achieve the best result in individual patients.

-These implants are used only as an aid to healing, they are not a substitute for normal intact tissue or bone, these implants are NOT intended to bear the weight of the body and early full weight bearing prior to complete bone healing may cause breakage of implant.

**IF EXCESSIVE LOADING CANNOT BE PREVENTED, AN IMPLANT SHOULD NOT BE USED.**

-The main goal of surgery with this implant is to establish bony fusion. Abnormal or excessive forces could lead to delayed union, non-union, or failure of the implant. Abnormal force loading and subsequent wear may be caused by:

- Uncorrected instability
- Improperly sized implant
- Inadequate soft tissue support
- Implant malposition
- Excessive motion
- Uncorrected or recurrent deformity
- Patient misuse or overactivity

-Proper fixation at the time of surgery is critical to the success of the procedure. Bone stock must be adequate to support the device.

-Over time, metallic implants may loosen, fracture, or cause pain after bone fracture or osteotomy is healed. Removal of metallic implants is at the surgeon's discretion, and the appropriateness of the selected procedure will be based on the surgeon's personal medical training and experience. It is imperative that adequate post-operative care and protection be provided by the surgeon.

-The implants are temporary internal fixation devices. The Internal fixation implants are designed to stabilize the Fracture bones during the normal healing process. After the fracture bone is fused, these devices serve no functional purpose and may be removed. Recovery time depends on the extent and type of fracture ( it may take from 4 to 12 months). Estimated functional life time of the devices is Approximately 4 to 12 months. After the bone heals, implants may be either left in place or removed, the final decision on implant removal is up to the surgeon.

#### **F.ADVVERSE EFFECTS**

In any surgical procedure, the potential for complications exists. The risks and complications with these implants include:

- 1- Infection or painful, swollen or inflamed implant site Fracture of the implant
- 2- Fracture,Loosening or dislocation of the implant requiring revision surgery
- 3- Allergic reaction(s) to implant material(s)
- 4- Migration of particle wear debris possibly resulting in a bodily response
- 5) Embolis

#### **G.HANDLING AND STERILIZATION**

This product is for single use only. Reuse of these devices may potentially result in serious patient harm. Examples of hazards related to the reuse of these devices include, but are not limited to: significant degradation in device performance, cross-infection, and contamination.

The Implants provided non-sterile must be sterilized by a hospital validated steam autoclaving process in appropriate protective wrapping when necessary. If necessary Components must be cleaned prior to sterilization in compliance with hospital validated cleaning process or Cleaning equipment manufacturers' user instructions and recommendations for chemical detergent is required.

The following process parameters are validated by Tipsan and recommended for sterilization.

**Steam Sterilization**

Cycle Type	Parameter	Minimum Set Point
Prevacuum 273°F (134°C)	Exposure Temperature	273°F (134°C)
	Exposure Time	5 minutes

These recommendations have been developed and validated using specific equipment. Due to variations in environment and equipment, it must be demonstrated that these recommendations produce sterility in your environment.

#### **H.STORAGE CONDITIONS**

All implants must be stored in a clean, dry environment and be protected from sunlight and extremes in temperature.

#### **I.MAGNETIC RESONANCE IMAGING (MRI) SAFETY**

There are inherent risks associated with the use of metallic implants in the MR environment; including component migration, heat induction, and signal interference or distortion near the component(s). Heat induction of metallic implants is a risk related to component geometry and material, as well as the MR power, duration, and pulse sequence. Since MR equipment is not standardized, the severity and likelihood of occurrence are unknown for these implants.

Tipsan implants have not been evaluated for safety and compatibility in the MR environment. Tipsan implants have not been tested for heating or migration in the MR environment. Since these devices have not been tested, Tipsan cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy.

These components are passive metallic devices, and as with all passive devices, there is potential for reciprocal interference with certain imaging modalities; including image distortion for MR and X-ray scatter in CT.

#### **J.SUBSTANCES CONTAINED BY THE DEVICE**

the devices do not contain or incorporate:

- a medicinal substance, including a human blood or plasma derivative, or
- tissues or cells, or their derivatives, of human origin, or
- tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No 722/2012; or
- substances which are carcinogenic, mutagenic or toxic to reproduction listed in Annex I, Chapter II, Section 10.4 of MDR 2017/745 .

#### **K.REMOVAL AND DISPOSAL**

Should a implant removal be necessary, the appropriate instrument should be used. Removed implants should be handled as bio-hazard products and be discarded in accordance with all local and national regulations for infection prevention and control.

#### **L.IMPLANT CARD**

As per MDR, Article 18 (3), Plates & Screws are exempted from the Implant Card obligation. A physical implant card is not delivered with Device.

#### **M.REPORTING OF SERIOUS INCIDENTS**

Surgeons/Health Facilities must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the serious incident has occurred;

#### **N. OTHER INFORMATION**

Access to “Summary of safety and clinical performance Report” :

Tipsan has drawn up a “Summary of Safety and Clinical Performance” report in compliance with Article 32 of 2017/745 Medical Device Regulation. The Summary of Safety and Clinical Performance (SSCP) can be accessed on the Eudamed public website (<https://ec.europa.eu/tools/eudamed>) or can be requested directly from Tipsan

TIPSAN TIBBİ ALETLER SAN. VE TİC. A.Ş.  
Kemalpaşa Cad.7404/1 Sk.No:3 Pınarbaşı- İzmir- Turkey  
Tel: +90 232 479 5654 www.tipsan.com.tr

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- (AB) 722/2012 sayılı Tüzük'te belirtildiği şekilde, hayvan kaynaklı dokular veya hücreler ya da onların türevlerini;