TIPSAN SPINAL PEDICULAR SCREW FIXATION SYSTEM

TPS-IFU-07

Rev.00, Issue Date:10.10.2022

The following languages are included in this packet:

The following languages are included in this packet:

English (EN) Turkish (TR)

For additional information and translations please contact the manufacturer or local distributor.

Attention Operating Surgeon
IMPORTANT MEDICAL INFORMATION
SPINAL PEDICULAR SCREW FIXATION SYSTEM

- A- DEFINITIONS
- B- DESCRIPTION
- C- INTENDED PURPOSE & INDICATIONS
- D- CONTRAINDICATIONS
- E- GENERAL PRECAUTIONS & WARNINGS
- F- ADVERSE EFFECTS
- G- HANDLING AND STERILIZATION
- H- STORAGE CONDITIONS
- I- MAGNATIC RESONANCE IMAGING (MRI) SAFETY
- J- SUBSTANCES CONTAINED BY THE DEVICE
- K - REMOVEL AND DISPOSAL
- L- IMPLANT CARD
- M- REPORTING OF SERIOUS INCIDENTS
- N- OTHER INFORMATION

A.DEFINITIONS:

Symbols and abbreviations may be used on the package label. The following table provides the definition of these symbols and abbreviations.

Table 1. Definitions of Symbols and Abbreviations

Symbol	Definition	Symbol	Definition
	Batch number		Date of manufacture
	Catalog number		Manufacturer
	Do not re-use		Unique Device Identifier Number
	Caution, consult accompanying documents		Medical Device
	Consult operating instructions		Non-sterile

B.DESCRPTION:

Tipsan Spinal Pedicle Screw System is a spinal internal fixation. It consists of various components including pedicle Screws, vertebral rods, connectors, Hooks, Occipital Plate.

Tipsan Spinal Pedicle Screw System includes following devices:

- Thoracic-Lumbar Spinal Pedicle Screw Fixation System
- Cervical Spinal Pedicle Screw Fixation System

Do not mix components from different manufacturer. All implantable devices are designed for single use only.

Materials:

The implants are manufactured from titanium alloy which conform to ASTM or ISO standards, or internal standards. The component material is provided on the package label

Table-2 : Component List & Materials & Application Method.

Components	Materials
All Components (Pedicle Screws, Rods, Connectors, Hooks, Occipital Plate)	ISO 5832-3 /Titanium Alloy

Implantation of the devices will expose patients to the following materials:

Material Description	Elements in Material	% Composition
Titanium alloy - Ti-6Al-4V (ISO 5832-3)	Aluminium	5,5 to 6,75
	Vanadium	3,5 to 4,5
	Iron	0,3 max.
	Oxygen	0,2 max.
	Carbon	0,08 max.

	Nitrogen	0,05 max.
	Hydrogen	0,015 max.
	Titanium	Balance

C. INTENDED PURPOSE & INDICATIONS

Intended Purpose:

-The Thoracic-Lumbar Spinal Pedicle Screw Fixation System is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the Thoracic, Lumbar Vertebra and Sacrum.

-Cervical Spinal Pedicle Screw Fixation System is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the cervical.

Indications:

Spinal Pedicle Screw Fixation System is indicated for

- 1) Degenerative disc disease (DDD)
- 2) Spondylolisthesis;
- 3) Trauma (i.e., fracture or dislocation);

Intended Clinical Performance:

Intended Clinical performance of device is:

- Pain reduction;
- spine stabilization and fusion

However, these performances are dependent on several parameters regardless of the device concerned:

- The patient (weight, activity level, expectations, age, etc.);
- The patient's bone quality
- The implantation quality.

Intended Patient populations:

The Device is designed to be used for skeletally mature patients. No limitations in regard to gender or aetiology

Intended user :

The Device is designated to be used only by Health care professional – Orthopaedic / Spinal Surgeon,

D. CONTRAINDICATIONS

Contraindications may be relative or absolute. The choice of a particular device must be carefully weighed against the patient's overall evaluation.

Circumstances listed below may reduce the chances of a successful outcome:

- 1) Any abnormality present which affects the normal process of bone remodeling including, but not limited to, severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis.
- 2) Insufficient quality or quantity of bone which would inhibit rigid device fixation.
- 3) Previous history of infection.
- 4) Excessive local inflammation.
- 5) Open wounds.
- 6) Any neuromuscular deficit which places an unsafe load level on the device during the healing period.
- 7) Obesity. An overweight or obese patient can produce loads on the spinal system which can lead to failure of the fixation of the device or to failure of the device itself. Obesity is defined according to the W.H.O. standards.
- 8) Patients having inadequate tissue coverage of the operative site.
- 9) A condition of senility, mental illness, or substance abuse. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications.
- 10) Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made prior to material selection or implantation.
- 11) Other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors, congenital abnormalities. These contraindications can be relative or absolute and must be taken into account by the physician when making his decision.

E. GENERAL PRECAUTIONS & WARNINGS

- The implantation of pedicle screw spinal systems must be performed only by experienced spinal surgeons having undergone the necessary specific training in the use of such systems because this is a technically demanding procedure presenting a risk of serious injury to the patient.

-The manufacturer is not liable for complications that might occur for reasons such as incorrect indications or surgical technique, unsuitable choice or handling of material and/or surgical instruments.

-The correct selection and sizing of the implant is extremely important. Selection of the proper size, shape, and design of the implant increases the potential for success. The implants require careful seating and adequate bone support.

-It is the responsibility of each surgeon using implants to consider the clinical and medical status of each patient and to be knowledgeable about all aspects of implant procedure and the potential complications that may occur. The benefits derived from implant surgery may not meet the patient's expectations or may deteriorate with time, necessitating revision surgery to replace the implant or to carry out alternative procedures. Revision surgeries with implants are common. The patient's mental status must also be considered. Willingness and/or ability to follow post-operative instructions may also impact the surgical outcome. Surgeons must balance many considerations to achieve the best result in individual patients.

The implants are temporary internal fixation devices.. After the Bone fusion, these devices serve no functional purpose and may be removed. Recovery time may take **from 4 to 12 months**. Estimated functional life time of the devices is Approximately 4 to 12 months. After the bone fusion, implants may be either left in place or removed, the final decision on implant removal is up to the surgeon.

- The choice of proper shape, size and design of the implant for each patient is crucial to the success of the surgery. The surgeon is responsible for this choice which depends on each patient. Patients who are overweight may be responsible for additional stresses and strains on the device which can speed up metal fatigue and/or lead to deformation or failure of the implants. The size and shape of the bone structures determine the size, shape and type of the implants. Once implanted, the implants are subjected to stresses and strains. These repeated stresses on the implants should be taken into consideration by the surgeon at the time of the choice of the implant, during implantation as well as in the post-operative follow-up period. Indeed, the stresses and strains on the implants may cause metal fatigue or fracture or deformation of the implants, before the bone graft has become completely consolidated. This may result in further side effects or necessitate the early removal of the osteosynthesis device.

- Improper selection, placement, positioning and fixation of these devices may result in unusual stress conditions reducing the service life of the implant. Contouring or bending of rods is recommended only if necessary according to the surgical technique of each system. Rods should only be contoured with the proper contouring instruments. Incorrectly contoured rods, or rods which have been repeatedly or excessively contoured must not be implanted. The surgeon is to be thoroughly familiar with the surgical procedure, instruments and implant characteristics prior to performing surgery.

- The surgeon must discuss all physical and psychological limitations inherent to the use of the device with the patient. This includes the rehabilitation regimen, physical therapy, and wearing an appropriate orthosis as prescribed by the physician. Particular discussion should be directed to the issues of premature weight bearing, activity levels, and the necessity for periodic medical follow-up.

-The surgeon must warn the patient of the surgical risks and made aware of possible adverse effects. The surgeon must warn the patient that the device cannot and does not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged as a result of strenuous activity or trauma, and that the device may need to be replaced in the future. If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g.,substantial walking, running, lifting, or muscle strain) the surgeon must advice the patient that resultant forces can cause failure of the device. Patients who smoke have been shown to have an increased incidence of non-unions. Surgeons must advise patients of this fact and warn of the potential consequences. For diseased patients with degenerative disease, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. In such cases, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief.

F.ADVVERSE EFFECTS

In any surgical procedure, the potential for complications exists. The risks and complications with these implants include:

- 1- Infection or painful, swollen or inflamed implant site Fracture of the implant
- 2- Fracture,Loosening or dislocation of the implant requiring revision surgery
- 3- Allergic reaction(s) to implant material(s)
- 4- Migration of particle wear debris possibly resulting in a bodily response
- 5- Embolis

G.HANDLING AND STERILIZATION

This product is for single use only. Reuse of these devices may potentially result in serious patient harm. Examples of hazards related to the reuse of these devices include, but are not limited to: significant degradation in device performance, cross-infection, and contamination.

The Implants provided non-sterile must be sterilized by a hospital validated steam autoclaving process in appropriate protective wrapping when necessary. If necessary Components must be cleaned prior to sterilization in compliance with hospital validated cleaning process or Cleaning equipment manufacturers' user instructions and recommendations for chemical detergent is required.

The following process parameters are validated by Tipsan and recommended for sterilization.

Steam Sterilization		
Cycle Type	Parameter	Minimum Set Point
Prevacuum 273°F (134°C)	Exposure Temperature	273°F (134°C)
	Exposure Time	5 minutes

These recommendations have been developed and validated using specific equipment. Due to variations in environment and

equipment, it must be demonstrated that these recommendations produce sterility in your environment.

H.STORAGE CONDITIONS

All implants must be stored in a clean, dry environment and be protected from sunlight and extremes in temperature.

I.MAGNETIC RESONANCE IMAGING (MRI) SAFETY

There are inherent risks associated with the use of metallic implants in the MR environment; including component migration, heat induction, and signal interference or distortion near the component(s). Heat induction of metallic implants is a risk related to component geometry and material, as well as the MR power, duration, and pulse sequence. Since MR equipment is not standardized, the severity and likelihood of occurrence are unknown for these implants.

Tipsan implants have not been evaluated for safety and compatibility in the MR environment. Tipsan implants have not been tested for heating or migration in the MR environment. Since these devices have not been tested, Tipsan cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy.

These components are passive metallic devices, and as with all passive devices, there is potential for reciprocal interference with certain imaging modalities; including image distortion for MR and X-ray scatter in CT.

J.SUBSTANCES CONTAINED BY THE DEVICE

the devices do not contain or incorporate:

- a medicinal substance, including a human blood or plasma derivative, or
- tissues or cells, or their derivatives, of human origin, or
- tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No 722/2012; or
- substances which are carcinogenic, mutagenic or toxic to reproduction listed in Annex I, Chapter II, Section 10.4 of MDR 2017/745 .

K.REMOVAL AND DISPOSAL

Should a implant removal be necessary, the appropriate instrument should be used. Removed implants should be handled as bio-hazard products and be discarded in accordance with all local and national regulations for infection prevention and control.

L.IMPLANT CARD

A physical implant card is delivered with Device, Health facilities / Hospitals should fill in this implant card and give it to Patients.

M.REPORTING OF SERIOUS INCIDENTS

Surgeons/Health Facilities must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the serious incident has occurred;

N. OTHER INFORMATION

Access to “Summary of safety and clinical performance Report” :

Tipsan has drawn up a “Summary of Safety and Clinical Performance” report in compliance with Article 32 of 2017/745 Medical Device Regulation. The Summary of Safety and Clinical Performance (SSCP) can be accessed on the Eudamed public website (<https://ec.europa.eu/tools/eudamed>) or can be requested directly from Tipsan

TIPSAN TIBBİ ALETLER SAN. VE TİC. A.Ş.
Kemalpaşa Cad.7404/1 Sk.No:3 Pınarbaşı- İzmir- Turkey
Tel: +90 232 479 5654 www.tipsan.com.tr

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- (AB) 722/2012 sayılı Tüzük'te belirtildiği şekilde, hayvan kaynaklı dokular veya hücreler ya da onların türevlerini;