TIPSAN TIBBİ ALETLER SAN. VE TİC. A.Ş.
Kemalpaşa Cad.7404/1 Sk.No:3 Pınarbaşı- İzmir- Turkey
Tel: +90 232 479 5654 www.tipsan.com.tr

TIPSAN HIP SYSTEMS

TPS-IFU-01 HIP

Rev.01, Issue Date:20.05.2024

The following languages are included in this packet:

The following languages are included in this packet:

English (EN) Turkish (TR)

For additional information and translations please contact the manufacturer or local distributor.

Attention Operating Surgeon
IMPORTANT MEDICAL INFORMATION
HIP SYSTEM

- A- DEFINITIONS
- B- DESCRIPTION
- C- INTENDED PURPOSE & INDICATIONS
- D- CONTRAINDICATIONS
- E- GENERAL PRECAUTIONS & WARNINGS
- F- ADVERSE EFFECTS
- G- HANDLING AND STERILIZATION
- H- STORAGE CONDITIONS
- I- MAGNATIC RESONANCE IMAGING (MRI) SAFETY
- J- SUBSTANCES CONTAINED BY THE DEVICE
- K - REMOVEL AND DISPOSAL
- L- IMPLANT CARD
- M- REPORTING OF SERIOUS INCIDENTS
- N- OTHER INFORMATION

A.DEFINITIONS:

Symbols and abbreviations may be used on the package label. The following table provides the definition of these symbols and abbreviations.

Table 1. Definitions of Symbols and Abbreviations

Symbol	Definition	Symbol	Definition
	Batch code		Date of manufacture
	Catalog number		Manufacturer
	Do not re-use		Sterilized using radiation
	Caution, consult accompanying documents		Do not use if package is damaged
	Consult operating instructions		Unique Device Identifier Number
	Use by		Medical Device
	Keep dry		Double sterile barrier system
	Keep away from sunlight		Contains hazardous substances
	Do not resterilise		

B.DESCRPTION:

Tipsan has a variety of hip joint replacement prostheses. The components for these systems include an acetabular component, femoral head, bipolar Head, femoral stem. These components can be utilized in a variety of configurations to assemble the final construct. Only components from Tipsan should be used to prevent mismatch or misalignment of components.

The implants are single use only devices.

Materials:

The femoral, acetabular, and cement restrictor components are manufactured from a variety of materials which include cobalt-chromium-molybdenum alloy, titanium alloy, Stainless steel, ceramic and ultra high molecular weight polyethylene (UHMWPE) all of which conform to ASTM or ISO standards, or internal standards. The component material is provided on the outside carton label

Table-2 : Hip System

Device Description	Cemented/Cementless Use - Primary/ Revision Surgery
MODULAR STRAIGHT STEM CEMENTED	Cemented / Primary Surgery
MODULAR STRAIGHT STEM UNCEMENTED	Uncemented / Primary Surgery
MODULAR BIPOLAR CUP	Uncemented / Primary Surgery
ACETABULAR CUP CEMENTED	Cemented / Primary Surgery
PRESSFIT CUP SYSTEM UNCEMENTED	Uncemented / Primary Surgery
CTD MODULAR STRAIGHT STEM CEMENTED	Cemented / Primary Surgery
CTPH MODULAR STRAIGHT STEM UNCEMENTED	Uncemented / Primary Surgery
PRIZMA MODULAR STEM UNCEMENTED	Uncemented / Primary Surgery
PRIZMA MODULAR REVISION STEM UNCEMENTED	Uncemented / Revision Surgery
T2 HIP PROSTHESIS	Uncemented / Primary Surgery
T2 REVISION HIP PROSTHESIS	Uncemented / Revision Surgery

Implantation of the devices will expose patients to the following materials:

Material Description	Elements in Material	% Composition
Titanium alloy - Ti-6Al-4V (ISO 5832-3)	Aluminium	5,5 to 6,75
	Vanadium	3,5 to 4,5
	Iron	0,3 max.
	Oxygen	0,2 max.
	Carbon	0,08 max.
	Nitrogen	0,05 max.
	Hydrogen	0,015 max.
	Titanium	Balance

Material Description	Elements in Material	% Composition
Co.Cr Alloy (ISO 5832-12)	Chromium	26,0 to 30,0
	Molybdenum	5,0 to 7,0
	Iron	0,75 max.
	Manganese	1,0 max.
	Silicon	1,0 max.
	Carbon	0,15 to 0,35
	Nickel	1,0 max.
	Nitrogen	0,25 max.
	Cobalt	Balance

Material Description	Elements in Material	% Composition
Casting Co.Cr Alloy Alloy (ISO 5832-4)	Chromium	26,5 to 30,0
	Molybdenum	4,5 to 7,0
	Nickel	1,0 max.
	Iron	1,0 max.
	Carbon	0,35 max.
	Manganese	1,0 max.
	Silicon	1,0 max.
	Cobalt	Balance

Material Description	Elements in Material	% Composition
UHMWPE (ISO 5834-2) (ISO 5834-1)	Ultra high molecular weight polyethylene	%100

Material Description	Elements in Material	% Composition (% mass Fraction)
Ceramic -yttria-stabilized tetragonal zirconia (Y-TZP) ISO 13356	Amount of $ZrO_2+HfO_2+Y_2O_3$	Min. %99
	Yttrium oxide Y_2O_3	4.5 to 6
	Hafnium oxide HfO_2	5 max
	Alumina Matrix Al_2O_3	0.5 Max.

	Other oxides	0.5 Max.
--	--------------	----------

Material Description	Elements in Material	% Composition
Stainless Steel (ISO 5832-1)	Carbon	0,030 max.
	Silicon	0,75 max.
	Manganese	2,0 max.
	Phosphorus	0,025 max.
	Sulfur	0,010 max.
	Nitrogen	0,10 max.
	Chromium	17,0 to 19,0
	Molybdenum	2,25 to 3,00
	Nickel	13,0 to 15,0
	Copper	0,50 max.
	Cobalt	0,10 max.
Iron	Balance	

Components manufactured by Cobalt Alloy material contain the following substances defined as CMR 1B in a concentration above 0.1% weight by weight:

-Cobalt; CAS No. 7440-48-4; EC No. 231-158-0.

Some preclinical and clinical evidence suggests that Co CA (cobalt-containing alloys) in medical devices are unlikely carcinogenic or reproductive hazards to patients.

C.INTENDED PURPOSE & INDICATIONS

Intended Purpose:

Tipsan hip systems are intended for use in hip arthroplasty for reduction or relief of pain and/or improved hip function in skeletally mature patients.

Indications:

Femoral Stem System and Acetabular Cup System (Total Hip Joint Replacement) are indicated for the following conditions :

- 1-non-inflammatory degenerative joint disease such as osteoarthritis, avascular necrosis, ankylosis, protrusio acetabuli, and painful hip dysplasia;
- 2-inflammatory degenerative joint disease such as rheumatoid arthritis;
- 3-correction of functional deformity; and,
- 4-revision procedures where other treatments or devices have failed (for Revision implant)

Bipolar System and Unipolar System (hemiarthroplasty -Partial Hip Joint Replacement) is indicated for the following conditions:

- 1-Pathological fractures of the femoral neck;
- 2-Non-union of femoral neck fractures;
- 3-Aseptic necrosis of the femoral head and neck; and,
- 4-Primary pathology in the young involving the femoral head but with a non deformed acetabulum.

Intended Clinical Performance:

Intended Clinical performance of device is:

- Pain reduction;
- The restoration of the functioning of the joint and the improvement of the mobility compared to the preoperative state;

However, these performances are dependent on several parameters regardless of the device concerned:

- The patient (weight, activity level, expectations, age, etc.);
- The patient's bone quality and the initial bone deficit;
- The quality of the soft tissues
- The implantation quality.

Intended Patient populations:

Tipsan Hip System is designed to be used for skeletally mature patients. No limitations in regard to gender

Intended user :

Tipsan Hip Systems are designated to be used only by Health care professional – Orthopaedic Surgeon,

D.CONTRAINDICATIONS

Patients should be warned of these contraindications.

Contraindications include:

- 1-overt infection;
- 2-distant foci of infections (which may cause hematogenous spread to the implant site);
- 3-rapid disease progression as manifested by joint destruction or bone absorption apparent on roentgenogram;
- 4-skeletally immature patients (patient is less than 21 years of age at the time of surgery);
- 5-cases where there is inadequate neuromuscular status (e.g., prior paralysis, fusion and/or inadequate abductor strength), poor bone stock, poor skin coverage around the joint which would make the procedure unjustifiable;
- 6-neuropathic joints;
- 7-hepatitis or HIV infection;
- 8-neurological or musculoskeletal disease that may adversely affect gait or weight-bearing.

E.GENERAL PRECAUTIONS & WARNINGS

The manufacturer is not liable for complications that might occur for reasons such as incorrect indications or surgical technique, unsuitable choice or handling of material and/or surgical instruments.

Preoperative

- The surgeon must evaluate each situation individually based on the patient's clinical presentation in making any decisions regarding implant selection. The surgeon must be thoroughly familiar with the implant, instruments and surgical procedure prior to performing surgery. The surgeon should contact Tipsan for product-specific surgical techniques.
- Patient selection should consider the following factors which could lead to increased risk of failure and can be critical to the eventual success of the procedure: the patient's weight, activity level, and occupation. Implant longevity and stability may be affected by these variables. A heavy-weight patient can produce high loads on the prosthesis, which can lead to failure of the prosthesis. The surgeon must consider the ability and willingness of the patient to follow instructions and to control their weight and activity level. Patients with high activity levels, poor bone quality, or heavyweight patients may not be candidates for a narrower femoral implant. Any joint replacement system, including the implant/bone interface, cannot be expected to withstand activity levels and loads as would normal healthy bone and will not be as strong, reliable, or durable as a natural human joint. The patient should not have unrealistic functional expectations for occupations or activities that include substantial walking, running, lifting, or muscle strain.
- Additional conditions presenting increased risk of failure include:
 - 1-uncooperative patient or patient with mental or neurologic conditions which can affect patient's ability or willingness to follow instructions;
 - 2-marked bone loss, severe osteoporosis, or revision procedures for which an adequate fit of the prosthesis cannot be achieved;
 - 3-metabolic disorders that may impair bone formation;
 - 4-osteomalacia;
 - 5-conditions that could impair or impede healing (e.g., alcohol or drug abuse, decubitus ulcer, end-stage diabetes, severe protein deficiency and/or malnutrition);
 - 6-pre-existing conditions commonly considered with any surgery including bleeding disorders, long-term steroidal therapy, immunosuppressive therapy, or high dosage radiation therapy.
- The patient should be warned of surgical risks, and made aware of possible adverse effects. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, has a finite expected service life, and may need to be replaced at some time in the future. The patient should also be advised of other risks that the surgeon believes should be disclosed. The patient should be advised that any noise or unusual sensation should be reported to the surgeon as it may indicate implant malfunction.
- Patients should be warned that the longevity of the implant may depend on their weight and level of activity
- Patient sensitivities or allergies to the materials of the implants, particularly to metal ions are possible. The type of materials for each implant is provided on the box label. The attending physician should include possible risks concerning sensitivities to implant materials as part of the preoperative planning.

Intraoperative

- Specialized instruments are available and must be used to assure the accurate implantation of prosthetic components. Do not mix instruments from different manufacturers. While rare, breakage of instruments may occur especially with extensive use or excessive force. For this reason, instruments should be examined for wear or damage prior to surgery.
- Do not mix components from different manufacturers. Failure to comply may result in implant failure and revision surgery
- Inspect devices prior to use for damage during shipment or storage or any out-of-box defects that might increase the likelihood of fragmentation during a procedure.
- **Correct selection of the prosthesis is extremely important.** Joint prostheses require careful seating and adequate bone support. Surgeons are encouraged to use their best medical judgment when choosing the proper implant size regardless of the endosteal area of the bone. Proper implant selection must consider design, fixation, patient weight, age, bone quality, size, activity level, preoperative level of health, and also the surgeon's experience and familiarity with the device. Implant longevity and stability may be affected by these variables. Surgeons should inform the patient about these factors.
- X-ray templates are used to estimate the size of the product to be used. The anatomy of the patient ultimately determines the size of the product for an individual patient. The extent of bone preparation is determined intraoperatively by reaming and/or broaching starting at the smallest size and continuing until bleeding cancellous bone is reached. Trial prostheses should be used to evaluate the position of the final implant and the joint range of motion. The final size of the implant selected during surgery may differ from the size originally planned during preoperative assessment or the combination chosen during preliminary trialing.
- **Cemented Application.** Care is to be taken to ensure complete support of all components of the prosthesis embedded in bone cement to prevent stress concentrations that may lead to failure of the procedure. Complete cleaning including complete removal of bone chips, bone cement fragments, and metallic debris, prior to closure of the prosthetic site is critical to prevent accelerated wear of the articular surfaces of the prosthesis. The Distal Centralizers are intended for use as part of a cemented total hip arthroplasty.
- **Non-Cemented Application.** Adequate fixation at the time of surgery is critical to the success of the procedure. Uncemented femoral stems and acetabular shells must press fit into the host bone, which necessitates precise operative technique and the use of specified instruments. Bone stock must be adequate to support the device.
- **NEVER** combine these metals in NON-ARTICULATING contact surfaces:
 - Stainless steel Modular Head / cobalt chrome alloy Femoral Stem

- Do not attempt to seat the implant beyond the envelope of femoral bone preparation. Forcing to seat the implant beyond the prepared femoral bone may increase the chance of bone fracture. In some cases, a portion of the proximal body with or without coating may be visible above the proximal resection level.
- The smaller sized femoral implants are intended for patients with narrower intramedullary femoral canals. The geometry of these implants is reduced to accommodate the anatomy of the narrower intramedullary femoral canal, which also decreases the fatigue-strength and load-bearing characteristics of the implant.
- Modular Components (Femoral Head and Stems, Modular Necks and Proximal Body). Scratching of femoral heads, modular necks and proximal and distal stem tapers should be avoided. Repeated assembly and disassembly of these components could compromise the locking action of the taper joint. Prior to assembly, surgical debris must be cleaned from the interior of the female seat of the proximal body to ensure proper locking. Ensure components are firmly seated to prevent disassociation. The femoral head, neck taper of the femoral component, modular neck tapers, body taper, female seat of the proximal body must be clean and dry before assembly.
- Cobalt chrome femoral heads are designed for use with cobalt-chromium-molybdenum, titanium alloy femoral stems.
- Fixation screws should be fully seated to ensure stable fixation of the shell, and avoid interference with the liner component. Before implanting, be certain the selected shell and liner are compatible. Prior to seating the liner component into the shell component, surgical debris must be cleaned from the interior of the shell and the shell must be thoroughly dried. Debris and fluid may inhibit the liner from locking into the shell component. Failure to properly seat the liner into the shell can lead to disassociation of the liner from the shell.
- Ceramic femoral heads and acetabular liners should not be placed on scratched or previously assembled metal tapers as this may lead to a ceramic fracture.
- Fracture of ceramic components is a serious complication. Special care must be taken with ceramic devices, which must not be used if dropped, even in the absence of any apparent damage. Only use a plastic tip to introduce the ceramic devices. Patients should be advised to report unusual noise and/or sharp pain as both can be an indication of fracture. Decision to revise should not be postponed as ceramic fragments can cause severe damage to surrounding soft tissue and metal components. Revision outcomes after ceramic fractures can be compromised by the remaining ceramic debris present in the tissue even after careful debridement. Damage has been reported in polyethylene and metal components used in revisions after ceramic fractures. Surgeons are advised to carefully consider all available implant options on an individual basis. It must be noted that removal of all components including femoral stems and acetabular shells may not prevent accelerated wear due to ceramic debris in the tissue. Partial or complete synovectomy has been recommended by some authors.

Postoperative

- While expected life time of devices is difficult to estimate, Device life time is Finite. Functional device life time depends on variable factors such as patient weight, activity, patient bone quality therefore It not possible to state the exact functional life time. According to the literature review about the performance and security of Joint replacement Implant devices at long term follow-up, Estimated functional life time of the devices is Approximately 10 to 15 years. This estimated lifetime may change depends on variable factors such as patient weight, activity, patient bone quality and surgeon experience.
- The patient must be advised of the limitations of the reconstruction and the need for protection of the prosthesis from full weight bearing until adequate fixation and healing have occurred. The patient should be cautioned to limit activities and protect the replaced joint from unreasonable stresses and possible loosening, fracture and/or wear, and follow the instructions of the physician with respect to follow-up care and treatment. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult.
- Periodic, long-term follow-up is recommended to monitor the position and state of the prosthetic components, as well as the condition of the adjoining bone. Periodic post-operative x-rays are recommended for close comparison with early post-op conditions to detect long term evidence of changes in position, loosening, bending, or cracking of components.
- The patient should be advised to report any pain, decrease in range of motion, swelling, fever, squeaking, clicking, popping, grating, or grinding noises, and unusual incidences. Patient reports of squeaking, clicking, popping, grating, or grinding should be carefully evaluated as they may indicate position changes in the components which may compromise the durability of the implants.

F.ADVERSE EFFECTS for hip arthroplasty implants can include:

- 1-Osteolysis (progressive bone resorption). Osteolysis can be asymptomatic and therefore routine periodic radiographic examination is vital to prevent any serious future complication.
- 2-Particulates leading to increased wear rates necessitating early revision.
- 3-Allergic reactions to materials; metal sensitivity; or reactions to wear debris that may lead to histological reactions, pseudotumor and aseptic lymphocytic vasculitis-associated lesions (ALVAL).
- 4-Delayed wound healing; Deep wound infection (early or late) which may necessitate removal of the prosthesis. On rare occasions, arthrodesis of the involved joint or amputation of the limb may be required.
- 5-A sudden drop in blood pressure intra-operatively due to the use of bone cement;
- 6-Damage to blood vessels or hematoma;
- 7-Temporary or permanent nerve damage, peripheral neuropathies and subclinical nerve damage as possible result of surgical trauma resulting in pain or numbness of the affected limb;
- 8-Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
- 9-Fatigue fracture of the prosthetic component can occur as a result of trauma, strenuous activity, improper alignment, incomplete implant seating, duration of service, loss of fixation, non-union, or excessive weight;
- 10-Dislocation, migration and/or subluxation of prosthetic components from improper positioning, trauma, loss of fixation and/or muscle and fibrous tissue laxity;
- 11-Periarticular calcification or ossification, with or without impediment to joint mobility;
- 12-Trochanteric non-union due to inadequate reattachment and or early weight bearing;
- 13-Trochanteric avulsion as a result of excess muscular tension, early weight bearing, or inadvertent intraoperative weakening;
- 14-Inadequate range of motion due to improper selection or positioning of components, by femoral impingement, and periarticular calcification;
- 15-Femoral or acetabular perforation or fracture; femoral fracture while seating the device; femoral fracture by trauma or excessive loading, particularly in the presence of poor bone stock;
- 16-Undesirable shortening or lengthening of the limb;
- 17-Aggravated problems of the affected limb or contralateral extremity by leg length discrepancy, excess femoral medialization, or muscle deficiency;
- 18-Pain.

G.HANDLING AND STERILIZATION

Implants are sterilized by gamma radiation. The immediate package label should be consulted for specific method of sterilization. Irradiated implants have been exposed to a minimum 25 and a maximum 40 kiloGrays of gamma radiation.

The products are supplied as sterile and should be considered sterile unless the inner package is opened or damaged. If the inner package integrity is compromised, contact the manufacturer for instructions. Remove from package, using aseptic OR technique, only after the correct size has been determined and the operative site has been prepared for final implantation. Always handle the product

with powder-free gloves, and avoid contact with hard objects that may damage the product. This is particularly important in handling porous coated and HA coated prostheses. Do not allow porous surfaces or HA surfaces to come in contact with cloth or other fiber releasing materials.

Devices labelled for single-use only should never be reused. Reuse of these devices may potentially result in serious patient harm. Examples of hazards related to the reuse of these devices include, but are not limited to: significant degradation in device performance, cross-infection, and contamination.

A prosthesis should never be resterilised or reused after contact with body tissues or fluids, but rather should be discarded. Tipsan does not take any responsibility for the use of implants resterilised after contact with body tissues or fluids.

WARNINGS:

- All packaging materials MUST be removed from the implant prior to implantation.
- Inspect packages for punctures or other damage prior to surgery, If the sterile barrier has been broken, return the component to Tipsan

H.STORAGE CONDITIONS

All implants must be stored in a clean, dry environment and be protected from sunlight and extremes in temperature.

I.MAGNETIC RESONANCE IMAGING (MRI) SAFETY

There are inherent risks associated with the use of metallic implants in the MR environment; including component migration, heat induction, and signal interference or distortion near the component(s). Heat induction of metallic implants is a risk related to component geometry and material, as well as the MR power, duration, and pulse sequence. Since MR equipment is not standardized, the severity and likelihood of occurrence are unknown for these implants.

Tipsan Hip Systems have not been evaluated for safety and compatibility in the MR environment. Tipsan Hip Systems have not been tested for heating or migration in the MR environment. Since these devices have not been tested, Tipsan cannot make a recommendation for the use of MRIs with these implants, neither for safety considerations nor imaging accuracy.

These components are passive metallic devices, and as with all passive devices, there is potential for reciprocal interference with certain imaging modalities; including image distortion for MR and X-ray scatter in CT.

J.SUBSTANCES CONTAINED BY THE DEVICE

the devices do not contain or incorporate:

- a medicinal substance, including a human blood or plasma derivative, or
- tissues or cells, or their derivatives, of human origin, or
- tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No 722/2012; or

K.REMOVAL AND DISPOSAL

Tipsan Implants are intended for permanent implantation. Should a implant removal be necessary, the appropriate instrument should be used. Removed implants should be handled as bio-hazard products and be discarded in accordance with all local and national regulations for infection prevention and control.

L.IMPLANT CARD

A physical implant card is delivered with Device, Health facilities / Hospitals should fill in this implant card and give it to Patients.

M.REPORTING OF SERIOUS INCIDENTS

Surgeons/Health Facilities must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the serious incident has occurred;

N. OTHER INFORMATION

Access to “Summary of safety and clinical performance Report” :

Tipsan has drawn up a “Summary of Safety and Clinical Performance” report in compliance with Article 32 of 2017/745 Medical Device Regulation. The Summary of Safety and Clinical Performance (SSCP) can be accessed on the Eudamed public website (<https://ec.europa.eu/tools/eudamed>) or can be requested directly from Tipsan

TIPSAN TIBBİ ALETLER SAN. VE TİC. A.Ş.
Kemalpaşa Cad.7404/1 Sk.No:3 Pınarbaşı- İzmir- Turkey
Tel: +90 232 479 5654 www.tipsan.com.tr

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- (AB) 722/2012 sayılı Tüzük'te belirtildiği şekilde, hayvan kaynaklı dokular veya hücreler ya da onların türevlerini;