



T08 PRIMARY KNEE MINIMAL INVASIVE

*the key in
orthopaedic
sciences*

FEMORAL COMPONENT CRUCIATE RETAINING (CR)

- CoCrMo Alloyed
- Cemented & Cementless
- 6 size, Left and Right



FEMORAL COMPONENT POSTERIOR STABILIZED (PS)

- CoCrMo Alloyed
- Cemented & Cementless
- 6 size, Left and Right



UHMWPE HIGH CROSSLINK INSERT

FIXED INSERT CRUCIATE RETAINING (CR)



- UHMWPE

FIXED INSERT POSTERIOR STABILIZED (PS)



- 6 Size Options

MOBILE INSERT CRUCIATE RETAINING (CR)



- 5 Different Thickness Sizes are Available

MOBILE INSERT POSTERIOR STABILIZED (PS)



FIXED TIBIAL COMPONENT

- CoCrMo Alloyed
- Cemented & Cementless
- 6 Size



MOBILE TIBIAL COMPONENT

- CoCrMo Alloyed
- Cemented & Cementless
- 6 Size



EXTENSION STEM

- CoCrMo Alloyed
- Extension Stem Connection
25 mm & 50 mm



PATELLA COMPONENT

- UHMWPE
- 5 Sizes
- 6 Sizes

